



# EE512 – Applied Biomedical Signal Processing

## Neural network architecture

Clémentine AGUET  
CSEM Signal Processing Group



- Labs



# Lab – Instructions

- Submit report as **single PDF file**
- Recommended to work in groups of **3 students**
- You can prepare one single report for the group (name1\_name2\_name3\_lab\_NN.pdf)
- But every member must upload the file on Moodle
- Python code is given and provided as **Jupyter notebooks**
- This practical session is not focused on coding but on questions testing your understanding and interpretation of the results.
- **2 exercises** in this lab session on real-life biomedical problems

# Labs – Exercise 1


- **ECG rhythm classification**

- GOAL: Train NN to classify different cardiac rhythms from single-lead ECG signals
- Data: <https://physionet.org/content/ecg-arrhythmia/1.0.0/>
- Cardiac rhythms:
  - Atrial fibrillation
  - Atrial flutter
  - Normal sinus rhythm
  - Sinus bradycardia
  - Sinus tachycardia
- 1500 single-lead (lead II) ECG signals of each rhythm

# Labs – Exercise 1

- ECG rhythm classification

Different in HR compared to normal sinus rhythm?  
HR irregularity?



# Labs – Exercise 1

- **ECG rhythm classification**

- Split data into training, validation and testing subsets stratified by rhythms
- 5 folds
  - 3 in training
  - 1 in validation
  - 1 in testing

| Subset     | Total | Atrial<br>fibrillation | Atrial flutter | Normal<br>sinus | Sinus<br>bradycardia | Sinus<br>tachycardia |
|------------|-------|------------------------|----------------|-----------------|----------------------|----------------------|
| Training   | 4500  | 900                    | 900            | 900             | 900                  | 900                  |
| Validation | 1500  | 300                    | 300            | 300             | 300                  | 300                  |
| Test       | 1500  | 300                    | 300            | 300             | 300                  | 300                  |

# Labs – Exercise 1

- **ECG rhythm classification**

- Scale ECG signals to have approximately unit variance
- Encode rhythm labels with one-hot encoding

|                     | <b>Atrial<br/>fibrillation</b> | <b>Atrial flutter</b> | <b>Normal<br/>sinus</b> | <b>Sinus<br/>bradycardia</b> | <b>Sinus<br/>tachycardia</b> |
|---------------------|--------------------------------|-----------------------|-------------------------|------------------------------|------------------------------|
| Atrial_fibrillation | 1.0                            | 0.0                   | 0.0                     | 0.0                          | 0.0                          |
| Sinus_bradycardia   | 0.0                            | 0.0                   | 0.0                     | 1.0                          | 0.0                          |
| Sinus_bradycardia   | 0.0                            | 0.0                   | 0.0                     | 1.0                          | 0.0                          |
| Atrial_flutter      | 0.0                            | 1.0                   | 0.0                     | 0.0                          | 0.0                          |
| Atrial_fibrillation | 1.0                            | 0.0                   | 0.0                     | 0.0                          | 0.0                          |
| Normal_sinus        | 0.0                            | 0.0                   | 1.0                     | 0.0                          | 0.0                          |
| Sinus_bradycardia   | 0.0                            | 0.0                   | 0.0                     | 0.0                          | 1.0                          |
| Atrial_flutter      | 0.0                            | 1.0                   | 0.0                     | 0.0                          | 0.0                          |

# Labs – Exercise 1

- **ECG rhythm classification**
  - CNN

```
class CnnModel(torch.nn.Module):  
  
    def __init__(self, input_shape, output_shape, kernel_size=5):  
        super().__init__()  
        self.input_shape = input_shape  
        self.output_shape = output_shape  
        self.layers = torch.nn.Sequential(  
            torch.nn.Conv1d(self.input_shape[0], 8, kernel_size, padding='same'),  
            torch.nn.BatchNorm1d(8),  
            torch.nn.ReLU(),  
            torch.nn.MaxPool1d(2),  
  
            torch.nn.Conv1d(8, 16, kernel_size, padding='same'),  
            torch.nn.BatchNorm1d(16),  
            torch.nn.ReLU(),  
            torch.nn.MaxPool1d(2),  
  
            torch.nn.Conv1d(16, 32, kernel_size, padding='same'),  
            torch.nn.BatchNorm1d(32),  
            torch.nn.ReLU(),  
            torch.nn.MaxPool1d(2),  
  
            torch.nn.Conv1d(32, 64, kernel_size, padding='same'),  
            torch.nn.BatchNorm1d(64),  
            torch.nn.ReLU(),  
            torch.nn.AdaptiveAvgPool1d(1),  
  
            torch.nn.Flatten(),  
            torch.nn.Linear(64, self.output_shape[0]),  
        )  
  
    def forward(self, x):  
        return self.layers(x)
```



# Labs – Exercise 1

- **ECG rhythm classification**

- CNN





|        | Name                                   | Type              | Params | In sizes     | Out sizes    |
|--------|--|-------------------|--------|--------------|--------------|
| 0      | model                                  | CnnModel          | 14.2 K | [1, 1, 1280] | [1, 5]       |
| 1      | model.layers                           | Sequential        | 14.2 K | [1, 1, 1280] | [1, 5]       |
| 2      | model.layers.0                         | Conv1d            | 48     | [1, 1, 1280] | [1, 8, 1280] |
| 3      | model.layers.1                         | BatchNorm1d       | 16     | [1, 8, 1280] | [1, 8, 1280] |
| 4      | model.layers.2                         | ReLU              | 0      | [1, 8, 1280] | [1, 8, 1280] |
| 5      | model.layers.3                         | MaxPool1d         | 0      | [1, 8, 1280] | [1, 8, 640]  |
| 6      | model.layers.4                         | Conv1d            | 656    | [1, 8, 640]  | [1, 16, 640] |
| 7      | model.layers.5                         | BatchNorm1d       | 32     | [1, 16, 640] | [1, 16, 640] |
| 8      | model.layers.6                         | ReLU              | 0      | [1, 16, 640] | [1, 16, 640] |
| 9      | model.layers.7                         | MaxPool1d         | 0      | [1, 16, 640] | [1, 16, 320] |
| 10     | model.layers.8                         | Conv1d            | 2.6 K  | [1, 16, 320] | [1, 32, 320] |
| 11     | model.layers.9                         | BatchNorm1d       | 64     | [1, 32, 320] | [1, 32, 320] |
| 12     | model.layers.10                        | ReLU              | 0      | [1, 32, 320] | [1, 32, 320] |
| 13     | model.layers.11                        | MaxPool1d         | 0      | [1, 32, 320] | [1, 32, 160] |
| 14     | model.layers.12                        | Conv1d            | 10.3 K | [1, 32, 160] | [1, 64, 160] |
| 15     | model.layers.13                        | BatchNorm1d       | 128    | [1, 64, 160] | [1, 64, 160] |
| 16     | model.layers.14                        | ReLU              | 0      | [1, 64, 160] | [1, 64, 160] |
| 17     | model.layers.15                        | AdaptiveAvgPool1d | 0      | [1, 64, 160] | [1, 64, 1]   |
| 18     | model.layers.16                        | Flatten           | 0      | [1, 64, 1]   | [1, 64]      |
| 19     | model.layers.17                        | Linear            | 325    | [1, 64]      | [1, 5]       |
| 14.2 K | Trainable params                       |                   |        |              |              |
| 0      | Non-trainable params                   |                   |        |              |              |
| 14.2 K | Total params                           |                   |        |              |              |
| 0.057  | Total estimated model params size (MB) |                   |        |              |              |

# Labs – Exercise 1

- **ECG rhythm classification**
  - Tensorboard

Overfitting?

Benefit for training longer?



# Labs – Exercise 1

- **ECG rhythm classification**
  - Evaluation with confusion matrices

Train set (accuracy = 0.764)

| References \ Predictions | atrial_fibrillation | atrial_flutter | normal_sinus_rhythm | sinus_bradycardia | sinus_tachycardia |
|--------------------------|---------------------|----------------|---------------------|-------------------|-------------------|
| atrial_fibrillation      | 609                 | 241            | 3                   | 36                | 11                |
| atrial_flutter           | 322                 | 503            | 4                   | 34                | 37                |
| normal_sinus_rhythm      | 28                  | 9              | 665                 | 142               | 56                |
| sinus_bradycardia        | 11                  | 6              | 30                  | 853               | 0                 |
| sinus_tachycardia        | 12                  | 40             | 36                  | 4                 | 808               |

Val set (accuracy = 0.732)

| References \ Predictions | atrial_fibrillation | atrial_flutter | normal_sinus_rhythm | sinus_bradycardia | sinus_tachycardia |
|--------------------------|---------------------|----------------|---------------------|-------------------|-------------------|
| atrial_fibrillation      | 194                 | 78             | 4                   | 21                | 3                 |
| atrial_flutter           | 124                 | 142            | 2                   | 28                | 4                 |
| normal_sinus_rhythm      | 13                  | 3              | 218                 | 45                | 21                |
| sinus_bradycardia        | 5                   | 2              | 15                  | 278               | 0                 |
| sinus_tachycardia        | 2                   | 18             | 12                  | 2                 | 266               |

Test set (accuracy = 0.716)

| References \ Predictions | atrial_fibrillation | atrial_flutter | normal_sinus_rhythm | sinus_bradycardia | sinus_tachycardia |
|--------------------------|---------------------|----------------|---------------------|-------------------|-------------------|
| atrial_fibrillation      | 191                 | 84             | 0                   | 20                | 5                 |
| atrial_flutter           | 172                 | 110            | 1                   | 12                | 5                 |
| normal_sinus_rhythm      | 7                   | 4              | 218                 | 48                | 23                |
| sinus_bradycardia        | 2                   | 0              | 9                   | 289               | 0                 |
| sinus_tachycardia        | 3                   | 16             | 14                  | 1                 | 266               |

Are some classes more difficult to classify?

# Labs – Exercise 1

- **ECG rhythm classification**
  - 2 custom architectures
  - Add layers: convolutional, pooling, FC, etc.

```
class CustomModel1(torch.nn.Module):  
  
    def __init__(self, input_shape, output_shape):  
        super().__init__()  
        self.input_shape = input_shape  
        self.output_shape = output_shape  
  
        # Implement you own model here.  
        self.layers = torch.nn.Sequential(  
            torch.nn.Flatten(),  
            torch.nn.Linear(np.prod(self.input_shape), self.output_shape[0]),  
        )  
  
    def forward(self, x):  
        return self.layers(x)
```

# Labs – Exercise 2

- **Heart rate classification**

- GOAL: Estimate heart rate (HR) from PPG and acceleration signal
- Data: <https://archive.ics.uci.edu/ml/datasets/PPG-DaLiA>
  - PPG signals
  - Acceleration signals
  - Reference HR computed from ECG signal
- Data collected during various activities but focus on:
  - Sitting
  - Walking



# Labs – Exercise 2

- **Heart rate classification**

- Preprocessing (already done)
  - Band-pass filtering between 0.4 and 4.0 Hz (24 – 240 bpm)
  - Resampling to 25 Hz
- Split signal into windows with overlap
  - Window length: 8 s
  - Shift length: 2 s



# Labs – Exercise 2

- Heart rate classification



# Labs – Exercise 2

- Heart rate classification





# Labs – Exercise 2

- **Heart rate classification**

- Extract signal windows with overlaps
  - Window length: 8 s
  - Shift length: 2 s
- 7420 windows with 1 or 4 channels
- Each window includes 200 samples (8 seconds at 25 Hz)
  - Input PPG (7420, 1, 200)
  - Input PPG + accelero (7420, 4, 200)

# Labs – Exercise 2

- **Heart rate classification**

- Split data into training, validation and testing subsets by subjects
  - 9 subjects in training set
  - 3 subjects in validation set
  - 3 subject in testing set

```
Subject used for training   : ['S1' 'S2' 'S3' 'S4' 'S5' 'S6' 'S7' 'S8' 'S9']  
Subject used for validation : ['S10' 'S11' 'S12']  
Subject used for testing    : ['S13' 'S14' 'S15']
```

18

- Scale windows to have approximately unit variance

# Labs – Exercise 2

- Heart rate classification

- CNN
- Input: PPG
- MSE loss function

PPG CNN config

```
{'model': {'input_shape': (1, 200),  
  'output_shape': (1,),  
  'n_convolutional_layers': 4,  
  'kernel_size': 5,  
  'n_initial_channels': 16,  
  'use_normalization': False,  
  'n_dense_layers': 3,  
  'n_units': 128,  
  'dropout': 0.0},  
  'optimizer': {'lr': 0.0001}}
```

|        | Name                                   | Type              | Params | In sizes     | Out sizes    |
|--------|--|-------------------|--------|--------------|--------------|
| 0      | model                                  | CnnModel          | 87.2 K | [1, 1, 200]  | [1, 1]       |
| 1      | model.layers                           | Sequential        | 87.2 K | [1, 1, 200]  | [1, 1]       |
| 2      | model.layers.0                         | Conv1d            | 96     | [1, 1, 200]  | [1, 16, 200] |
| 3      | model.layers.1                         | ReLU              | 0      | [1, 16, 200] | [1, 16, 200] |
| 4      | model.layers.2                         | MaxPool1d         | 0      | [1, 16, 200] | [1, 16, 100] |
| 5      | model.layers.3                         | Conv1d            | 2.6 K  | [1, 16, 100] | [1, 32, 100] |
| 6      | model.layers.4                         | ReLU              | 0      | [1, 32, 100] | [1, 32, 100] |
| 7      | model.layers.5                         | MaxPool1d         | 0      | [1, 32, 100] | [1, 32, 50]  |
| 8      | model.layers.6                         | Conv1d            | 10.3 K | [1, 32, 50]  | [1, 64, 50]  |
| 9      | model.layers.7                         | ReLU              | 0      | [1, 64, 50]  | [1, 64, 50]  |
| 10     | model.layers.8                         | MaxPool1d         | 0      | [1, 64, 50]  | [1, 64, 25]  |
| 11     | model.layers.9                         | Conv1d            | 41.1 K | [1, 64, 25]  | [1, 128, 25] |
| 12     | model.layers.10                        | ReLU              | 0      | [1, 128, 25] | [1, 128, 25] |
| 13     | model.layers.11                        | AdaptiveAvgPool1d | 0      | [1, 128, 25] | [1, 128, 1]  |
| 14     | model.layers.12                        | Flatten           | 0      | [1, 128, 1]  | [1, 128]     |
| 15     | model.layers.13                        | Linear            | 16.5 K | [1, 128]     | [1, 128]     |
| 16     | model.layers.14                        | ReLU              | 0      | [1, 128]     | [1, 128]     |
| 17     | model.layers.15                        | Linear            | 16.5 K | [1, 128]     | [1, 128]     |
| 18     | model.layers.16                        | ReLU              | 0      | [1, 128]     | [1, 128]     |
| 19     | model.layers.17                        | Linear            | 129    | [1, 128]     | [1, 1]       |
| 87.2 K | Trainable params                       |                   |        |              |              |
| 0      | Non-trainable params                   |                   |        |              |              |
| 87.2 K | Total params                           |                   |        |              |              |
| 0.349  | Total estimated model params size (MB) |                   |        |              |              |

# Labs – Exercise 2

- **Heart rate classification**

- CNN
- Input: PPG + accelero
- MSE loss function

PPG ACC CNN config



```
{'model': {'input_shape': (4, 200),
  'output_shape': (1,),
  'n_convolutional_layers': 4,
  'kernel_size': 5,
  'n_initial_channels': 16,
  'use_normalization': False,
  'n_dense_layers': 3,
  'n_units': 128,
  'dropout': 0.0},
  'optimizer': {'lr': 0.0001}}
```

|        | Name                                   | Type              | Params | In sizes     | Out sizes    |
|--------|--|-------------------|--------|--------------|--------------|
| 0      | model                                  | CnnModel          | 87.5 K | [1, 4, 200]  | [1, 1]       |
| 1      | model.layers                           | Sequential        | 87.5 K | [1, 4, 200]  | [1, 1]       |
| 2      | model.layers.0                         | Conv1d            | 336    | [1, 4, 200]  | [1, 16, 200] |
| 3      | model.layers.1                         | ReLU              | 0      | [1, 16, 200] | [1, 16, 200] |
| 4      | model.layers.2                         | MaxPool1d         | 0      | [1, 16, 200] | [1, 16, 100] |
| 5      | model.layers.3                         | Conv1d            | 2.6 K  | [1, 16, 100] | [1, 32, 100] |
| 6      | model.layers.4                         | ReLU              | 0      | [1, 32, 100] | [1, 32, 100] |
| 7      | model.layers.5                         | MaxPool1d         | 0      | [1, 32, 100] | [1, 32, 50]  |
| 8      | model.layers.6                         | Conv1d            | 10.3 K | [1, 32, 50]  | [1, 64, 50]  |
| 9      | model.layers.7                         | ReLU              | 0      | [1, 64, 50]  | [1, 64, 50]  |
| 10     | model.layers.8                         | MaxPool1d         | 0      | [1, 64, 50]  | [1, 64, 25]  |
| 11     | model.layers.9                         | Conv1d            | 41.1 K | [1, 64, 25]  | [1, 128, 25] |
| 12     | model.layers.10                        | ReLU              | 0      | [1, 128, 25] | [1, 128, 25] |
| 13     | model.layers.11                        | AdaptiveAvgPool1d | 0      | [1, 128, 25] | [1, 128, 1]  |
| 14     | model.layers.12                        | Flatten           | 0      | [1, 128, 1]  | [1, 128]     |
| 15     | model.layers.13                        | Linear            | 16.5 K | [1, 128]     | [1, 128]     |
| 16     | model.layers.14                        | ReLU              | 0      | [1, 128]     | [1, 128]     |
| 17     | model.layers.15                        | Linear            | 16.5 K | [1, 128]     | [1, 128]     |
| 18     | model.layers.16                        | ReLU              | 0      | [1, 128]     | [1, 128]     |
| 19     | model.layers.17                        | Linear            | 129    | [1, 128]     | [1, 1]       |
| 87.5 K | Trainable params                       |                   |        |              |              |
| 0      | Non-trainable params                   |                   |        |              |              |
| 87.5 K | Total params                           |                   |        |              |              |
| 0.350  | Total estimated model params size (MB) |                   |        |              |              |

# Labs – Exercise 2

- **Heart rate classification**

- HR prediction for validation set (with or without accelero)



Difference between with or without acceleration signals?

# Labs – Exercise 2

- Heart rate classification

- CNN
- Input: PPG + accelero
- MSE loss function
- Add batch normalization after each convolution layer

```
PPG ACC norm CNN config
{'model': {'input_shape': (4, 200),
  'output_shape': (1,),
  'n_convolutional_layers': 4,
  'kernel_size': 5,
  'n_initial_channels': 16,
  'use_normalization': True,
  'n_dense_layers': 3,
  'n_units': 128,
  'dropout': 0.0},
 'optimizer': {'lr': 0.0001}}
```



|        | Name                                   | Type              | Params | In sizes     | Out sizes    |
|--------|--|-------------------|--------|--------------|--------------|
| 0      | model                                  | CnnModel          | 88.0 K | [1, 4, 200]  | [1, 1]       |
| 1      | model.layers                           | Sequential        | 88.0 K | [1, 4, 200]  | [1, 1]       |
| 2      | model.layers.0                         | Conv1d            | 336    | [1, 4, 200]  | [1, 16, 200] |
| 3      | model.layers.1                         | BatchNorm1d       | 32     | [1, 16, 200] | [1, 16, 200] |
| 4      | model.layers.2                         | ReLU              | 0      | [1, 16, 200] | [1, 16, 200] |
| 5      | model.layers.3                         | MaxPool1d         | 0      | [1, 16, 200] | [1, 16, 100] |
| 6      | model.layers.4                         | Conv1d            | 2.6 K  | [1, 16, 100] | [1, 32, 100] |
| 7      | model.layers.5                         | BatchNorm1d       | 64     | [1, 32, 100] | [1, 32, 100] |
| 8      | model.layers.6                         | ReLU              | 0      | [1, 32, 100] | [1, 32, 100] |
| 9      | model.layers.7                         | MaxPool1d         | 0      | [1, 32, 100] | [1, 32, 50]  |
| 10     | model.layers.8                         | Conv1d            | 10.3 K | [1, 32, 50]  | [1, 64, 50]  |
| 11     | model.layers.9                         | BatchNorm1d       | 128    | [1, 64, 50]  | [1, 64, 50]  |
| 12     | model.layers.10                        | ReLU              | 0      | [1, 64, 50]  | [1, 64, 50]  |
| 13     | model.layers.11                        | MaxPool1d         | 0      | [1, 64, 50]  | [1, 64, 25]  |
| 14     | model.layers.12                        | Conv1d            | 41.1 K | [1, 64, 25]  | [1, 128, 25] |
| 15     | model.layers.13                        | BatchNorm1d       | 256    | [1, 128, 25] | [1, 128, 25] |
| 16     | model.layers.14                        | ReLU              | 0      | [1, 128, 25] | [1, 128, 25] |
| 17     | model.layers.15                        | AdaptiveAvgPool1d | 0      | [1, 128, 25] | [1, 128, 1]  |
| 18     | model.layers.16                        | Flatten           | 0      | [1, 128, 1]  | [1, 128]     |
| 19     | model.layers.17                        | Linear            | 16.5 K | [1, 128]     | [1, 128]     |
| 20     | model.layers.18                        | ReLU              | 0      | [1, 128]     | [1, 128]     |
| 21     | model.layers.19                        | Linear            | 16.5 K | [1, 128]     | [1, 128]     |
| 22     | model.layers.20                        | ReLU              | 0      | [1, 128]     | [1, 128]     |
| 23     | model.layers.21                        | Linear            | 129    | [1, 128]     | [1, 1]       |
|        |  |                   |        |              |              |
| 88.0 K | Trainable params                       |                   |        |              |              |
| 0      | Non-trainable params                   |                   |        |              |              |
| 88.0 K | Total params                           |                   |        |              |              |
| 0.352  | Total estimated model params size (MB) |                   |        |              |              |



# Labs – Exercise 2

- **Heart rate classification**

- HR prediction for validation set (with or without batch normalization)



Difference between with or without batch normalization?

# Labs – Exercise 2

- **Heart rate classification**

- CNN
- Input: PPG + accelero
- MSE loss function
- Batch normalization after each convolution layer
- Add dropout after each dense layer (except last one)

```
PPG ACC norm dropout CNN config
{'model': {'input_shape': (4, 200),
  'output_shape': (1,)},
  'n_convolutional_layers': 4,
  'kernel_size': 5,
  'n_initial_channels': 16,
  'use_normalization': True,
  'n_dense_layers': 3,
  'n_units': 128,
  'dropout': 0.5},
  'optimizer': {'lr': 0.0001}}
```



|        | Name                                   | Type              | Params | In sizes     | Out sizes    |
|--------|--|-------------------|--------|--------------|--------------|
| 0      | model                                  | CnnModel          | 88.0 K | [1, 4, 200]  | [1, 1]       |
| 1      | model.layers                           | Sequential        | 88.0 K | [1, 4, 200]  | [1, 1]       |
| 2      | model.layers.0                         | Conv1d            | 336    | [1, 4, 200]  | [1, 16, 200] |
| 3      | model.layers.1                         | BatchNorm1d       | 32     | [1, 16, 200] | [1, 16, 200] |
| 4      | model.layers.2                         | ReLU              | 0      | [1, 16, 200] | [1, 16, 200] |
| 5      | model.layers.3                         | MaxPool1d         | 0      | [1, 16, 200] | [1, 16, 100] |
| 6      | model.layers.4                         | Conv1d            | 2.6 K  | [1, 16, 100] | [1, 32, 100] |
| 7      | model.layers.5                         | BatchNorm1d       | 64     | [1, 32, 100] | [1, 32, 100] |
| 8      | model.layers.6                         | ReLU              | 0      | [1, 32, 100] | [1, 32, 100] |
| 9      | model.layers.7                         | MaxPool1d         | 0      | [1, 32, 100] | [1, 32, 50]  |
| 10     | model.layers.8                         | Conv1d            | 10.3 K | [1, 32, 50]  | [1, 64, 50]  |
| 11     | model.layers.9                         | BatchNorm1d       | 128    | [1, 64, 50]  | [1, 64, 50]  |
| 12     | model.layers.10                        | ReLU              | 0      | [1, 64, 50]  | [1, 64, 50]  |
| 13     | model.layers.11                        | MaxPool1d         | 0      | [1, 64, 50]  | [1, 64, 25]  |
| 14     | model.layers.12                        | Conv1d            | 41.1 K | [1, 64, 25]  | [1, 128, 25] |
| 15     | model.layers.13                        | BatchNorm1d       | 256    | [1, 128, 25] | [1, 128, 25] |
| 16     | model.layers.14                        | ReLU              | 0      | [1, 128, 25] | [1, 128, 25] |
| 17     | model.layers.15                        | AdaptiveAvgPool1d | 0      | [1, 128, 25] | [1, 128, 1]  |
| 18     | model.layers.16                        | Flatten           | 0      | [1, 128, 1]  | [1, 128]     |
| 19     | model.layers.17                        | Linear            | 16.5 K | [1, 128]     | [1, 128]     |
| 20     | model.layers.18                        | ReLU              | 0      | [1, 128]     | [1, 128]     |
| 21     | model.layers.19                        | Dropout           | 0      | [1, 128]     | [1, 128]     |
| 22     | model.layers.20                        | Linear            | 16.5 K | [1, 128]     | [1, 128]     |
| 23     | model.layers.21                        | ReLU              | 0      | [1, 128]     | [1, 128]     |
| 24     | model.layers.22                        | Dropout           | 0      | [1, 128]     | [1, 128]     |
| 25     | model.layers.23                        | Linear            | 129    | [1, 128]     | [1, 1]       |
|        |  |                   |        |              |              |
| 88.0 K | Trainable params                       |                   |        |              |              |
| 0      | Non-trainable params                   |                   |        |              |              |
| 88.0 K | Total params                           |                   |        |              |              |
| 0.352  | Total estimated model params size (MB) |                   |        |              |              |



# Labs – Exercise 2

- **Heart rate classification**

- HR prediction for validation set (with or without dropout)



# Labs – Exercise 2

- **Heart rate classification**

- Tensorboard

ppg\_cnn\version\_0







ppg\_acc\_cnn\version\_0



ppg\_acc\_norm\_cnn\version\_0



ppg\_acc\_norm\_dropout\_cnn\version\_0



# Labs – Exercise 2

- **Heart rate classification**

- Evaluation

| MSE | Model                            | Training | Validation | Test   |
|-----|----------------------------------|----------|------------|--------|
|     | PPG                              | 253.15   | 518.95     | 199.11 |
|     | PPG + ACC                        | 127.88   | 263.46     | 198.26 |
|     | PPG + ACC + batch norm           | 26.12    | 232.19     | 71.65  |
|     | PPG + ACC + batch norm + dropout | 39.89    | 239.53     | 49.74  |
| MAE | Model                            | Training | Validation | Test   |
|     | PPG                              | 12.37    | 19.01      | 12.43  |
|     | PPG + ACC                        | 8.78     | 13.10      | 12.05  |
|     | PPG + ACC + batch norm           | 3.75     | 11.48      | 6.16   |
|     | PPG + ACC + batch norm + dropout | 4.68     | 11.72      | 5.07   |